



Item Name : Adapters Socket To Cone Tee Connection

Item Code : 040.2008.01

Item Description :

Adapters Socket To Cone Tee Connection – Laboratory Glassware

Adapters Socket To Cone Tee Connection are precision-engineered laboratory connectors designed for controlled fluid distribution across applications.

Key Features

- **Tee Connection Design** – Allows simultaneous fluid transfer between three glassware components.
- **Socket to Cone Adaptability** – Ensures compatibility with various laboratory setups.
- **Leak-Proof & Contamination Resistant** – Maintains sample integrity during laboratory procedures.
- **Standardized Joint Compatibility** – Available in multiple socket and cone sizes for seamless integration.


Why Choose Our Adapters Socket To Cone Tee Connection?

Designed for laboratory precision, these adapters enhance fluid transfer efficiency while maintaining airtight sealing. Their tee configuration

Applications

- **Chemical & Pharmaceutical Research** – Supports controlled multi-port fluid transfer in experimental setups.
- **Biotech & Medical Laboratories** – Ensures efficient sample containment for analytical procedures.
- **Industrial & Scientific Glassware Integration** – Facilitates fluid transfer in manufacturing and testing environments.

Item Image :



Cat. No.	Socket	Cone
040.2008.01	14/23	14/23

040.2008.02	14/23	19/26
040.2008.03	14/23	24/29
040.2008.04	14/23	29/32
040.2008.05	19/26	19/26
040.2008.06	19/26	24/29
040.2008.07	19/26	29/32
040.2008.08	24/29	24/29
040.2008.09	24/29	29/32
040.2008.10	29/32	29/32

44, Vikas puri, P.O.Ind.Estate, Jagadhri Road,
A/Cantt-133 006. Haryana, India.

+0171 – 2698107

info@glacierindia.com
<https://www.glacierindia.com>